












COLOR METHOD	4 - 19
Color Method Before & After	6
All-In-One Kits	7
Build Your Own Tub	8
Tub Accessories	9
Trays Covers	10
Multi-Mod Racks	11
Cassettes & Instrument Holders	12 - 13
Steri-Containers	14
Bur & Endo Guards/Blocks	15 - 17
Rings & Tape	18
Consulting Services	19
SOLATION 2	0 - 25
Mr. Thirsty® One-Step	20 - 21
Airway Armor®	22
Pink Petal [™]	23
Insti-Dam [®]	24
Lingua-Fix® SE Cushion	25
MOUTH MIRRORS	26 - 31
Crystal HD® Mirrors	. 26 - 29
Mirror Magic [®] Anti-Fog Solution	30 - 31
OPERATORY/MISC 3	32 - 42
Denture & Retainer Boxes	32 - 33
Drawer Organizers	34 - 35
Steri-Soaker and CAD/CAM Block Organizer	36
E-Z Shelf & Accessories	37
Cabinet Trays & Lab Pans	38
HVE & SE Accessories and Dispensers	39
E-Z Storage and Handi-Hopper	40
Cuspidor/Dispos-A-Screen Traps and Needle Capper	41
Mighty Mixer Bowls, Spatula, and Organizers	42





Small Changes, Big Results

Zirc's Color Method creates a streamlined, stress-free day-to-day routine from setup to sterilization. Implementing this simple system adds a visual element to both instrument and material management yielding a workflow properly optimized for efficiency, leaving your team more confident and less stressed.



BEFORE **CLUTTERED DRAWERS CAUSES STRESS AND CHAOS**







STAY EFFICIENT WITH COLOR-CODED ORGANIZATION.



All-In-One Kits

_

Zirc's procedure tubs are designed to store and organize materials to keep your inventory under control.



Syringe Composite Tub



Capsule Composite Tub



Endo Tub



Crown Prep Tub



Hygiene Tub



Oral Surgery Tub



COMPLETE TUB

(20Z455) Available in all colors

Includes:

- Procedure Tub
- Tub Liner
- Tub Dividers
- Single Tub Cup/Cover

- Double Tub Cup/Cover
- Divided
 Slide Tray
- Safe-Lok®
 Cover for Tub

Dimensions: 32.07 x 28.57 x 10.16cm



Composite Kits s

Store and organize composite materials with our convenient Composite Kits.



SYRINGE COMPOSITE KIT (20Z489)

Syringe holder allows you to store up to 40 syringe tubes.

Includes:

- Syringe 20-Unit Stand
- 1 Long Tub Cup/Cover
- 1 Divided Slide Tray
- 1 Syringe Clip-On Holder
- Shade Label Sheet 2pk



CAPSULE COMPOSITE KIT (20Z479)

Long Tub Cups provide a total of 20 compartments for capsule storage.

Includes:

- Tub Liner
- 4 Long Tub Cup/Cover
- 1 Divided Slide Tray
- 1 Capsule Clip-On Holder
- Shade Label Sheet 2pk

Build Your Own Tub



PROCEDURE TUB (20Z463) Available in all colors

Allows you to "build" your own setup.

Dimensions: 31.11cm x 27.62cm x 6.98cm





SAFE-LOK® TUB COVER (20Z459)

Locking cover for securing materials. Keeps materials protected during storage and transportation.

Dimensions: 32.07cm x 28.57cm x 3.65cm





TUB LINER (20Z462)

Foam liner that sits in bottom of Tub. Tub not included. Replace as needed.

Dimensions: 24.45cm x 22.54cm x .32cm



TUB BOX (20**Z474)** Gray only

Designed to hold/organize supplies. Non-slip bottom for secure placement. Flip-top cover.

Dimensions: 9.52cm x 9.21cm x 4.76cm





TUB DIVIDERS (20Z461) Gray only

Creates up to 24 compartments. May be cut to create larger compartments.

Dimensions: 24.45cm x 21.91cm x 3.33cm







TUB INSERT (20Z465) Gray only

10 compartments for easy organization. Sits on rim of Tub, allowing for storage underneath.

Dimensions: 27.3cm x 24.61cm x 2.54cm



* Products not included

Tub Cups with Covers | s

Fill Tub Cups with brushes, clamps, cotton pellets, prophy angle tips, prophy paste, pumice, syringe tips, implant tools, etc. Designed to be used inside the Procedure Tub.



LONG TUB CUP WITH COVER (20Z473) Gray only

Internal dividers are removable.

Dimensions: 21.27cm x 5.87cm x 2.86cm



DOUBLE TUB CUP WITH COVER (20Z472) Gray only

Dimensions: 8.57cm x 5.71cm x 2.86cm



SINGLE TUB CUP WITH COVER (20Z471) Gray only

Dimensions: 5.71cm x 4.13cm x 2.86cm



Tub Accessories

Sits on rim of Tub and allows for storage underneath.



DIVIDED SLIDE TRAY (20Z466) Clear

Dimensions: 24.76cm x 0.95cm x 2.7cm



FLAT SLIDE TRAY (20Z467) Clear

Dimensions: 24.76cm x 0.95cm x 2.7cm



Preprinted labels for marking composite materials within the Tub organizers.

Dimensions: 17.78cm x 10.25cm x 0.04cm



SYRINGE 20-UNIT STAND (20Z487)

Stands up for easy access to syringe-able materials. Lays flat to fit inside Tub or drawer. Holds up to 20 syringes. Has space for 2 Syringe 10-Unit Add-Ons. (Stand can hold up to 40 syringes.)

Dimensions: 22.22cm x 12.1cm x 3.97cm





SYRINGE CLIP-ON HOLDER (20Z486)

Allows quick access to composite during procedure. Conveniently clips on any side of the Tub or B-Lok Tray. Holds 15 syringes and 3 capsules.

Dimensions: 20.32cm x 6.19cm x 1.9cm







SYRINGE 10-UNIT ADD-ON (20Z488)

An Add-On for the Syringe 20-Unit Stand. Gives additional space for syringe storage. Holds 10 syringes.

Dimensions: 21.75cm x 5.87cm x 1.90cm





CAPSULE CLIP-ON HOLDER (20Z476)

Allows quick access to composite during procedure. Conveniently clips on any side of the Tub or B-Lok Tray. Holds 8 capsules, and 4 syringes.

Dimensions: 20.32cm x 5.08cm x 1.43cm



Our exclusive B-Lok trays and covers provide an excellent method of infection control and safety that allows for two benefits other trays do not:

- Prevents any contaminated items from falling off of the used tray during transportation from operatory to central sterilization. The contaminated items stay enclosed until you are ready for the cleaning and sterilization process.
- Set up your procedure trays ahead of time with all instrumentation needed. The Safe-Lok® cover securely locks to the tray, not only keeping items clean, but also allowing for safe transportation to chairside. The cover is removed once the patient is in the chair and the procedure begins.





B-LOK FLAT TRAY (20Z401) Available in all colors

Wide edges for secure grip. Non-slip feet for stability. High-gloss surface helps prevent staining.

Dimensions: 33.97cm x 24.45cm x 2.22cm







Holds 12 hand instruments. Bur holder for 9 FG, RA, or HP burs. Comes with Insti-Grip. Non-slip feet for stability. Self-lubricating, high-gloss surface helps prevent staining. Holder for Syringe Tip or Saliva Ejector.

Dimensions: 33.97cm x 24.45cm x 2.22cm





MINI TRAY (20Z101) Available in all colors

Wide edges for secure grip. Non-slip feet for stability. High-gloss surface helps prevent staining. Mini Trays are great for hygiene/exam setup.

Dimensions: 23.81cm x 16.19cm x 2.22cm





INSTI-GRIP (20Z454I)

To be used with the B-Lok divided Tray only. Holds instruments quietly and safely in place.

Dimensions: 16.83cm x 0.63cm x 0.95cm

for b-lok divided tray



I I

SAFE-LOK® TRAY COVER **B-SIZE**

(20Z445)

Fits B-Lok Flat Tray and Divided Trays only. Slide locking mechanism securely locks cover to tray.

Dimensions: 35.24cm x 25.08cm x 1.90cm



SAFE-LOK® TRAY COVER **TALL B-SIZE**

(20Z446)

The Tall size accommodates our Performance Tool Organizer and almost any stainless steel cassette. Fits B-Lok Flat and Divided Trays only. Side locking mechanism securely locks cover to tray.

Dimensions: 35.24cm x 25.08cm x 3.49cm



SAFE-LOK® MINI TRAY COVER F-SIZE

(207103)

Side locking mechanism securely locks cover to tray.

Dimensions: 25.08cm x 16.83cm x 2.54cm





TRAY COVER B-SIZE (NON-LOCKING) (20Z441)

Dimensions: 35.56cm x 25.72cm x 0.79cm



MINI TRAY COVER F-SIZE (NON-LOCKING)

Dimensions: 24.45cm x 16.83cm x 2.22cm









The Safe-Lok® Tray Cover helps in compliance with OSHA and CDC, providing safe transportation back to the central sterilization area.

Are you in compliance?

OSHA I Immediately or as soon as possible after use, contaminated reusable sharp shall be placed in a puncture-resistant, leak-proof and labeled or color-coded containers until properly reprocessed. (To avoid spillage of contents, it is suggested that they be covered and secured prior to moving.)

CDC Contaminated instruments must be "contained" prior to transporting them. Containment can be defined as instrument cassettes or cages in which dental care items are cleaned, sterilized, and stored until point of use. Sharp instruments should not be carried openly to the instrument-processing area. Percutaneous injuries may occur to other dental healthcare personnel (DHCP) or patients in hallways leading to the instrument-processing area.

Multi-Mod Racks Adjustable to fit your needs!

Tub and tray racks are effective additions to a dental practice. These racks have non-slip feet so they are stable on any counter or cabinet surface. With a chrome-plated finish, these racks can fit any dental practice design.

Adjusts to hold Mini Trays

MULTI-MOD RACK 8-PLACE (217410)

Holds 8 Trays or 4 Tubs with covers. Sturdy metal construction; chrome-plated. Non-slip feet. Adjustable to hold trays and tubs.

S





MULTI-MOD RACK 6-PLACE (21Z105)

Holds 6 trays or 3 tubs with covers. Adjustable to hold trays and tubs. Sturdy metal construction; chrome-plated. Non-slip feet.

S

Why Use Zirc's Cassettes?

A low-cost investment with a quick return! Individual instrument handling takes time away from treating patients and increases overall reprocessing time. Zirc's cassettes improves instrument efficiency, safety, and compliance without breaking the bank.

Product Benefits

- Zirc's cassettes will save you time and help prevent injury.
- · Keep instrument setups together; no sorting.
- Protects your instruments from damage and prolongs their life.

191mm x 330mm 150mm

- · Prevents the team from a potential exposure.
- Better patient perception, making patients less anxious.

Product Features

- · Designed for superior water flow and cleaning open top, round holes, and radius corners.
- Keeps instrument setups together through sterilization and procedure processes.
- Non-slip feet for stability and safety.
- AUTOCLAVABLE.
- Available in 16 colors for color coding by procedure.

Sterilization Pouch Size

E-Z Jett® 5-Place

E-Z Jett® 8-Place	191mm x 330mm	200mm
E-Z Jett® 8-Place Slim	191mm x 330mm	200mm
E-Z Jett [®] 10-Place	254mm x 356mm	250mm
Performance Tool Organizer	254mm x 356mm	250mm
-175		
		-
	15	
	1	
		2/3
		13
		1
		333
	1	- 7
		3
	15	
		- 277
		-










Cassettes & Instrument Holders

_



E-Z JETT[®] CASSETTE 5-PLACE

(50Z925) Available in all colors

Holds 5 instruments. Great for basic exam setup! A/W Syringe Tip Holder included. Covers stay attached to the cassette and fold back easily. Fits Statim Sterilizer.

Dimensions: 18.57cm x 7.94cm x 3.02cm





E-Z JETT[®] CASSETTE 10-PLACE

(50Z930) Available in all colors

Holds 10 instruments. Accessory compartment in the middle for extra storage. Covers stay attached to the cassette and fold back easily.

Dimensions: 20.16cm x 17.78cm x 2.96cm





E-Z JETT° CASSETTE 8-PLACE

(50Z927) Available in all colors

Holds 8 instruments. Covers stay attached to the cassette and fold back easily. A/W Syringe Tip Holder included

Dimensions: 20.16cm x 11.27cm x 2.86cm





INSTRUMENT MAT (20Z464) Available in all colors

Holds 10 hand instruments. Reverse side holds 4 surgical instruments. Non-slip, soft material. Organizes instruments in order of procedure use. Can be cut to fit Mini Tray.

Dimensions: 17.14cm x 5.08cm x 0.95cm





Size comparison with the 8-Place E-Z Jett® Cassette




E-Z JETT° CASSETTE SLIM 8-PLACE

(50Z911) Available in all colors

Holds 8 instruments. Covers stay attached to the cassette and fold back easily. Compact size. Fits the Statim. A/W Syringe Tip Holder included.

Dimensions: 17.46cm x 9.05cm x 3.02cm





MINI MAT (20Z480) Available in all colors

Holds 8 instruments. Organizes instruments at chairside and for sterilization. Non-slip feet for stability.

Dimensions: 9.84cm x 1.27cm x 0.79cm

INSTRUMENT STRAP (50Z917) Gray only

Locking strap secures instruments.

Dimensions: 9.21cm x 1.27cm x 0.31cm





A deeper, roomier cassette made for large and specialty instruments. Designed for superior water flow for superb cleaning. Compact styling to fit in most sterilizers. Covers stay attached to the cassette and fold back easily. Comes with two strips of 10 customizable removable feet to design your own layout.

Ideal for:

- Surgical elevators and forceps
- Orthodontic pliers
- Hemostats
- Scissors
- · Any instrument that does not fit in your regular cassette

PERFORMANCE TOOL ORGANIZER

(50Z940) Available in all colors

Dimensions: 20.8cm x 10.79cm x 3.09cm

PERFORMANCE TOOL ORGANIZER **REPLACEMENT FEET**

(50Z944) Gray only

Two strips of 10.











Closed



STERI-CONTAINER

(50Z900) Available in all colors

Keeps instrument setups together through sterilization processes. Holds 10-14 instruments.

Dimensions: 20.64cm x 4.76cm x 4.76cm







COMPACT STERI-CONTAINER

(50Z905) Available in all colors

Easy open, push down latch. Keeps instrument setups together through sterilization processes. Statim 2000 holds 7 Compact Steri-Containers. Holds 8-12 instruments.

Dimensions: 18.10cm x 3.81cm x 3.81cm













Bur Blocks

Magnetic Bur Blocks 500 | 本」



- Fits FG, RA, and HP Burs
- Magnetic base holds burs securely
- · Non-slip feet for stability
- Drain holes help prevent rusting of burs



Recessed area to hold used burs.



MAGNETIC BUR BLOCK 8-HOLE

(50Z400) Available in all colors

Covers sold separately.

Dimensions: 3.81cm x 3.81cm x 1.27cm



MAGNETIC BUR BLOCK 14-HOLE

(50Z403) Available in all colors

Covers sold separately.

Dimensions: 7.46cm x 3.81cm x 1.27cm

Bur Block Covers s



- · Covers available to maintain cleanliness of burs during storage
- 2.54cm covers for FG and RA burs; 5.08cm cover (14-hole only) for lab and surgical burs
- Covers are NOT autoclavable
- Surface disinfect only



MAGNETIC BUR BLOCK 8-HOLE COVER 2.54CM HIGH STANDARD (50Z401)

Dimensions: 4.13cm x 4.13cm x 2.7cm



MAGNETIC BUR BLOCK 14-HOLE COVER 2.54CM HIGH STANDARD (50Z404)

Dimensions: 7.94cm x 4.13cm x 2.7cm



MAGNETIC BUR BLOCK 14-HOLE COVER 5.08CM HIGH HP LAB BURS (50Z405)

Dimensions: 7.94cm x 4.13cm x 5.71cm


WHITE **GREEN** BEIGE GRAY PLUM MIDNIGHT COPPER N. BLUE N. YELLOW N. GREEN N. ORANGE N. PURPLE

Steri-Bur Guards & Adapters

This gadget has it all! Not only does it protect your burs during cleaning and sterilization, but it can also be positioned at three different angles to give you easy access at chairside.



- Fits FG, RA, and HP Burs.
- Adjustable height for burs of different lengths (12-Hole Tall Bur Guard not adjustable).
- Holds burs in position for maximum cleaning effectiveness.
- Drain holes help prevent rusting of burs.






STERI-BUR GUARD 12-HOLE (50Z406) Available in all colors

Accommodates burs from 12mm-30mm.

Dimensions: 7.3cm x 0.95cm x 3.49cm

UNIVERSAL BUR ADAPTER 5PK (50Z412)

Can be used with 12- or 22-Hole Bur Guard. Enables long and short burs to be used in the same bur guard without short burs falling out.



a short bur adapter is included to stop any shorter burs from falling out



STERI-BUR GUARD TALL 12-HOLE (50Z408I) Gray only

Accommodates burs from 30.9mm-37.1mm.

Dimensions: 7.3cm x 0.95cm x 4.13cm



STERI-BUR GUARD 22-HOLE

(50Z410) Available in all colors

Accommodates burs from 12mm-30mm.

Dimensions: 12.38cm x 1.11cm x 3.81cm




Steri-Endo Guards & Holders

_

A unique delivery system to organize, store, and sterilize endodontic files. The Endo Guard swivels 360° with stops at every 45° for operator's choice. Large slotting provides complete visibility and draining during cleaning.

Steri-Endo Guard

- Holds most lengths for staging, cleaning, and sterilizing files and reamers.
- Millimeter gauge measures and positions rubber stopper.
- Counter on locking bar tracks the number of sterilizations.
- · Right- or left-hand use; optional angling.
- · Available in Hand or Engine styles.
- Single-use foam insert cleans/holds files (packaged with 12).






STERI-ENDO GUARD ENGINE (50Z452) Available in all colors

Organizes up to 11 engine and 5 hand files/reamers.

Dimensions: 13.65cm x 0.95cm x 5.56cm



STERI-ENDO GUARD HAND (50Z450) Available in all colors

Organizes up to 16 hand files/reamers.

Dimensions: 13.65cm x 0.95cm x 5.56cm

FOAM INSERTS 48PK (50Z451) Single use

FOAM INSERTS 144PK (50Z451V) Single use



ENDO ASSIST (50Z460) Available in all colors

Built-in millimeter gauge measures and positions rubber stopper. For rightor left-hand use. Small and large rings are included to accommodate most finger sizes. Flat bottom enables the Endo Assist to stand freely on its own. Single-use foam insert cleans/holds files (packaged with 12).

Dimensions: 3.81cm x 4.44cm x 4.13cm



ENDO ASSIST STAND (50Z462) Black only

Holds endo files and reamers in upright position. Foam provides fast and easy cleaning of the files during use.

Dimensions: 4.76cm x 4.76cm x 3.17cm





Dimensions: 26.67cm x 19.68cm x 2.22cm

ENDO ASSIST STAND FOAM 144PK (50Z466V) *Single use*

Dimensions: 26.67cm x 19.68cm x 2.22cm



Rings & Tape

How to effectively identify your instrument setup

1 color for procedure

Vibrant Orange is for identifying the procedure

1 color for practitioner (optional)

Vibrant Blue shows ownership to the practitioner



The diagonal setup of the orange rings shows the order of how the instruments should be used. An optional 2nd color may be used to identify the practitioner.




Place your instrument to determine which ring size you need.



E-Z ID Rings 💹 | 淅 |

- Prepackaged in an easy-to-use dispenser
- · Available in Small, Large, and Extra Large sizes to fit all instrument handles



25PK

Available in all colors

- (M) (70Z100)
- (G) (70Z200)
- (70Z000)

VIBRANT SYSTEM 200PK

Assorted colors (N, O, P, R, Q, and S)

- (M) (70Z110)
- (G) (70Z210)
- (70Z010)

JEWEL SYSTEM 200PK

Assorted colors (E, J, T, I, M, and U)

- (M) (70Z115)
- (G) (70Z215)
- (70Z015)

CLASSIC SYSTEM 200PK

Assorted colors (A, B, D, G, and I)

- (M) (70Z105)
- (G) (70Z205)
- (70Z005)



E-Z ID Tape [136]

- Tape dispenser has built-in tape cutter
- Dispenser roller prevents tape from twisting



E-Z ID TAPE ROLL (70Z300) Available in all colors

Includes 1 - 3.05m roll.

E-Z ID TAPE SYSTEM

Includes 8 - 0.91m roll.

Jewel (70Z315)

Assorted colors (E, J, T, I, M, and U)

Vibrant

(70Z310)

Assorted colors (N, O, P, Q, R, and S)

Classic

(70Z305)

Assorted colors (A, B, D, G, and I)







Consulting Services

_

Color Method will rescue your team from clinical clutter and disorganization. It's time to introduce a color-coded organization system to your practice. Schedule a consultation or RESET and we'll do it for you!

PLAN

1

GUIDANCE

Learn how Color Method will change the way you work! In this plan, you will receive a free one-on-one discussion with a Clinical Efficiency Consultant. Let's discuss your goals and challenges and get your questions answered!

2

DEEPER ASSESSMENT

We'll dive deeper into your clinical processes and send you a sample kit so you can experience the products for yourself. We'll also draw up a customized Color Method proposal including all the Zirc products you need and photos of how to organize your materials. Here's what's included:

- A Product Sample Kit
 - Inside is over 25 of Zirc's organizational products for you to explore
- A High-Level Assessment
- Our consultants will evaluate your clinical processes
- Final Assessment & Proposal Outlining
 - Current potential blindspots
 - Recommendations for improvement
- Your custom Color Method proposal

3

THE ULTIMATE RESET

RESET is a one-day workshop in your practice. Our Clinical Efficiency Consultants work side-by-side with your team to restructure your instruments and materials. We'll bring all the Zirc products with us and personally introduce you to the award-winning Color Method!



Mr. Thirsty® One-Step

Hands-free retraction, isolation, and evacuation provided by one simple device! Attach Mr. Thirsty® to your existing HVE, perform the procedure, and dispose when done! Rather than time spent on sterilizing or tracking number of uses, it can be spent on additional procedures, daily tasks, or catching a break! According to a time and motion study, Mr. Thirsty® can reduce chair procedure time by up to 29%!



mr. thirsty

- Built-in bite block for patient comfort
- Connects directly to HVE valve
- Keeps work environment dry
- Fully assembled
- Single use
- · Latex-free

new smaller size is perfect for smaller mouths & pediatric patients

VARIETY PACK

16mm Adaptor (50Z98816)

10mm Adaptor (50Z98810)

Includes 10 of each size:

- (10) Small
- (10) Medium
- (10) Large









With One Step, Achieve Simplicity



Affordable no equipment investment



Disposableno time wasted in
maintenance or sterilization



Adaptable connects to existing HVE, allowing it to be used in any operatory



Trimmable

available in three sizes that can be trimmed to fit, no need to overstock with a variety of sizes



Comfortable

provides tongue and cheek retraction and includes a built-in bite block, creating comfort for you and your patient



50PK

16mm Adaptor

- (50Z989SM16) Small
- (50Z989MD16) Medium
- (50Z989LG16) Large

10mm Adaptor

- (50Z989SM10) Small
- (50Z989MD10) Medium
- (50Z989LG10) Large

100PK

16mm Adaptor

- (50Z989SM-B16) Small
- (50Z989MD-B16) Medium
- (50Z989LG-B16) Large

10mm Adaptor

- (50Z989SM-B10) Small
- (50Z989MD-B10) Medium
- (50Z989LG-B10) Large

10MM/16MM ADAPTORS 1PK

10MM Adaptor (50Z993) 16mm Adaptor (50Z994)






With all the small items used in dentistry, (burs, files, implants and veneers) teams always run the risk of something being aspirated or ingested. Airway Armor® was created to eliminate the unnecessary risk for patients and protection for the practice.



3 convenient sizes



Ideal for Crown & Bridge, Endodontics, Extractions, Implants, Inlays/Onlays, Orthodontics, Veneers

- Aspiration and ingestion prevention
- Patient protection and comfort
- · Easy placement and removal
- 3 convenient sizes small, medium, and large
- Single use
- Latex-free







15PK TRIAL KIT (50Z960TRI)

Includes 5 of each size:

- **(5)** Small
- (5) Medium
- (5) Large

AIRWAY ARMOR® 25PK

(50Z960SM) Small (50Z960MD) Medium (50Z960LG) Large



Pink Petal[™]

Created by hygienist Karrie Velky, the Pink Petal[™] is a great product for all dental practices looking to fulfill the need of hands-free, low volume evacuation. Easily slide the Pink Petal™ on any saliva ejector, creating a device that will provide continuous suction.



A wonderful adjunct tool to free up a clinician's hand

- Holds low-volume suction tip in place
- Protects and retracts the cheek
- Patient can close around suction with ease
- · Easily slides onto saliva ejector
- Soft, comfortable material
- Single use
- Latex-free





Insti-Dam®

Features for BOTH Relaxed Fit and Latex-Free Insti-Dam®

- Built-in flexible medium gauge material that comes pre-assembled with frame and pre-punched hole
- Pre-punched hole makes positioning easy and allows eyes and nose to be open and comfortable
- · Additional holes can be punched
- · Compact design fits outside patients' lips, easy-to-use, non-threatening, and comfortable for patients
- · Radiographs may be taken by folding it to the side. No need to remove for mid procedure x-ray
- · Single use
- · Latex-free
- 44.8% faster placement than traditional dam.

INSTI-DAM®





LATEX-FREE INSTI-DAM® RELAXED FIT 20PK (50Z457)

Ideal for posterior isolation. The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing, or pulling.

Dimensions: 3.2mm thickness x 10.79cm diameter

Patent #7044738



LATEX-FREE INSTI-DAM® 20PK (50Z459)

Ideal for anterior isolation. Made with polyisoprene.

Dimensions: 3.2mm thickness x 10.79 diameter



INSTI-DAM® DISPENSER

(50Z471) Only available in white

Can be mounted to wall or cabinet or can sit on a counter. Non-slip rubber bottom for stability. Holds 35 Insti-Dams®. Made for Latex-Free Insti-Dam® ONLY.

Dimensions: 12.7cm x 13.33cm x 15.24cm

S





Lingua-Fix®

_

- - · Lingua-Fix® provides a dry work area
 - Holds tongue in a steady position to protect from injuries
 - Simple operation connects to low-volume suction valve
 - Comfortable, no sharp edges or corners
- Large suction surface avoids blockage and is more comfortable
- Reduces sedimentation of drilling byproducts in pipes and tubing

suction or evacuation

- Singe use
- · Comes pre-assembled

LINGUA-FIX® 10PK (50Z953S)

LINGUA-FIX® 50PK (50Z953)

SE Cushion

_

· Used with a Saliva Ejector

• Prevents tissue grabbing and tissue irritation

Increases patient comfort

• Quiets the noise of suction

• Single use



SE CUSHION 100PK (50Z970)

helpful tip ideal for use with rubber dams

WHITE BLUE GREEN BEIGE GRAY PLUM TEAL RED MIDNIGHT COPPER N. BLUE N. YELLOW N. GREEN N. ORANGE N. PURPLE N. PIN A B D G I E J M T U N O P Q R S

High Definition Imagery

What a view! The brighter the workspace, the better the illumination and visibility. With Zirc's Crystal HD® lens, you are visually working in a 40% brighter environment compared to the Rhodium lens that is currently the industry standard.

If you haven't made the change to HD - what are you waiting for?









Crystal HD® Mouth Mirrors

_


Why Use Zirc's Resin Mirrors?





⁻ Dr. Dimitris Tsanaktsidis, Microscopic & Aesthetic Dentistry

CUTTING EDGE IMAGERY



Color Reflectivity (In Wave Lengths)

Rhodium = Reflects color 75% accurately

Crystal HD® = Reflects color 99% accurately



Image provided by Dr. Lori Trost

Crystal HD® Mouth Mirrors

What's Your Mirror Size?





CRYSTAL HD® MOUTH MIRROR SOFT GRIP BLACK/GRAY



12PK 4PK

#4 | (50Z363) #4 | (50Z361) #5 | (50Z364) #5 | (50Z362)



CRYSTAL HD® MOUTH MIRROR SOFT GRIP 12PK



Assorted Vibrant						
2 of each color:						
N, O, P, Q, R, S						

E, J, T, I, M, U #4 | (50Z363-JWL) #4 | (50Z363-NEON)

Assorted Jewel

2 of each color:

#5 | (50Z364-NEON) #5 | (50Z364-JWL)



CRYSTAL HD® MOUTH MIRROR THIN GRIP 12PK



Assorted Vibrant **Assorted Jewel** 2 of each color: 2 of each color: N, O, P, Q, R, S E, J, T, I, M, U

#4 | (50Z353-NEON) #4 | (50Z353-JWL) #5 | (50Z354-NEON) #5 | (50Z354-JWL)



MEGA Mirrors

MEGA MIRROR

(50Z309) Color code: N, P, Q, R, S

A fun toy/prop for your patients. 6.35cm lens and 39.37cm in length.

S



Crystal HD® Mirror Heads & Handles

_



CRYSTAL HD° RESIN MOUTH MIRROR SCREW-ON HEADS (SIMPLE STEM) 12PK



One colorAssorted VibrantAvailable in all2 of each color:colors and K (black)N, O, P, Q, R, S

#4 | (50Z392) #4 | (50Z392-NEON) #5 | (50Z393) #5 | (50Z393-NEON)

Assorted Jewel 2 of each color: 2 of each color: E, J, T, I, M, U Assorted Classic 2 of each color: A, B, D, G, I, K

#4 | (50Z392-JWL) #4 | (50Z392-PSTL) #5 | (50Z393-JWL) #5 | (50Z393-PSTL)



CRYSTAL HD® STAINLESS MOUTH MIRROR SCREW-ON HEAD (SIMPLE STEM)



Stainless Steel 12pk #3 | (50R389) #4 | (50R390) #5 | (50R391)



Stainless Steel Head #0 | (50R388)





CRYSTAL HD® RESIN MOUTH MIRROR DOUBLE-SIDED HEADS (SIMPLE STEM)



Assorted Gray
Vibrant Head Mirror Head
4pk 4pk

1 of each color: N. O. R. S #4 | (50Z394)

#4 | (50Z394-NEON)

Crystal HD® Handles

_





ERGO GRIP HANDLE (SIMPLE STEM)

(50Z347) Available in K (black) & all colors except A

Dimensions: 12mm diameter x 13.33mm length

Mirror Magic[®] Anti-Fog System

The Mirror Magic[®] Anti-Fog System prevents fogging on mouth mirrors and offers the practitioner a quick method for removing debris during procedures. The disposable swipe pads stick to the back of a gloved hand for easy access. The Anti-Fog System requires only a small amount of solution that provides fog-free working time and is economic in size, making it easy to store.





The mirror magic® Anti-Fog System decreases fogging while also decreasing the build-up of debris on my mirror. This affords me outstanding visualization, which is crucial to performing quality dentistry. mirror magic® has truly made a difference in my practice and is one of those features I didn't realize I needed; and now I can't imagine practicing without it.

- Andrew M. Goldsmith, DDS, DICOI, FIALD

MIRROR MAGIC® TRIAL KIT (50R200X)

2 - 6ml bottles of Anti-Fog Solution. 24 Swipe Pads.



Without mirror magic®



With mirror magic®



THE PERFECT PAIR

The superior visual acuity and 40% brighter Crystal HD® lens combined with the quick debris removal and Anti-Fog Solution of mirror magic® have proven to be the perfect pair. Increase chairside efficiency and speed up your process with this dynamic duo.







MIRROR MAGIC* ANTI-FOG SOLUTION 10PK (50R205)

Each bottle contains 6ml each and lasts approx. 50 patients.



MIRROR MAGIC® SWIPE PADS 96PK (50R210) Single use



CLARITY COLLECTION (50R220)

Includes: Anti-Fog Solution 10pk, Swipe Pads 96pk, 1 Crystal HD* Mouth Mirror, 1 mirror magic* Organizer.





MIRROR MAGIC® ORGANIZER (50R201)

Holds: 48 Swipe Pads and 1-6ml bottle of Anti-Fog Solution. Store on operatory counter or in drawer. Non-slip bottom for secure placement.

Dimensions: 9.52cm x 9.21cm x 4.76cm

Note: box only, does not include mirror magic®

S



Denture & Retainer Boxes



Retainer & Denture Boxes S

One piece construction. Solid locking mechanism.

RETAINER BOX 12PK

2.54CM DEEP (25X500)

Dimensions: 7.94cm x 6.98cm x 2.54cm

3.81CM DEEP (25X550)

Dimensions: 7.94cm x 6.98cm x 3.81cm

Single Colors

TEAL	NEON BLUE	NEON GREEN	NEON PURPLE	NEON YELLOW	NEON PINK
	N	P	P	0	S

Assorted Pack

contains 3 of each color below





Dimensions: 11.11cm x 8.57cm x 5.08cm

Single

Assorted Pack contains 3 of each color below







Personalized Imprinting

Our **new and improved** imprinting service is quick and easy. Upload your own logo or select a full color icon from our pre-made Fun or Professional sets. Add up to four lines of text for your business name, phone number, website, or your own unique message. It's as easy as that!

IMPRINTED RETAINER BOXES 24PK

One piece construction. Solid locking mechanism.

2.54CM DEEP 3.81CM DEEP (25P500) (25P550)

Dimensions: 7.94cm x 6.98cm x 2.54cm
Dimensions: 7.94cm x 6.98cm x 3.81cm

IMPRINTED RETAINER BOXES 1.000PK

2.54CM DEEP (25PB500J) 3.81CM DEEP (25PB550J)

Dimensions: 7.94cm x 6.98cm x 2.54cm

Dimensions: 7.94cm x 6.98cm x 3.81cm

Single Colors



Assorted Pack

contains 6 of each for 24pk & 250 of each for 1000pk





IMPRINTED DENTURE BOX 24PK (30P800)

Dimensions: 11.11cm x 8.57cm x 5.08cm

IMPRINTED DENTURE BOX 1,000PK (30PB800I)

Dimensions: 11.11cm x 8.57cm x 5.08cm

Single Colors



Assorted Pack

contains 6 of each for 24pk & 250 of each for 1000pk



ORDER AT ZIRC.COM

WHITE BLUE GREEN BEIGE GRAY PLUM TEAL RED MIDNIGHT COPPER N. BLUE N. YELLOW N. GREEN N. ORANGE N. PURPLE N. PINK A B D G I E J M T U N O P Q R S

Drawer Organizer

Zirc is excited to announce a new product **coming soon** – a modular Drawer Organizer. All Zirc products are designed with time in mind for an organized and efficient dental practice, this new Drawer Organizer is no different.







DRAWER ORGANIZER

(20Z210) Color Code: A, G, I, N, P, R, S

Each Drawer Organizer will come with 1 large divider, 2 medium dividers and 4 small dividers. Additional dividers can be purchased.

Dimensions: 29.21cm x 16.76cm x 5.08cm

Large Divider (Included)

Dimensions: 28.84cm x 5cm x 2.42cm

DRAWER ORGANIZER DIVIDERS

Color Code: A, G, I, N, P, R, S

Dividers may be ordered separate. Size (S) comes in 2 pieces and size (M) in 1 piece

(20Z207) Small

Dimensions: 8.1cm x 5cm x 1.59cm

(20Z208) Medium

Dimensions: 16.42cm x 5cm x 1.59cm





















BEIGE G

The Steri-Soaker

An efficient, high-volume germicide tray designed for maximum use of chemicals to clean/disinfect instruments. It is designed to hold instruments for easy cleaning and is large enough to accommodate Zirc Steri-Containers, Cassettes, Bur Blocks, and Bur Guards.





STERI-SOAKER (50Z805) Color code: B, I

Calendar on lid conveniently reminds staff when to change the solution; use color code rings to mark dates. Transfer basket's ribbed and perforated design allows quick drain-off and easy pick-up of instruments. Holds 3.78L of liquid; fill lines shown in quarts and liters. Hinged fitted cover.

Dimensions: 35.24cm x 24.45cm x 13.02cm

Safe to use chemicals: cleaners (enzymes, detergents), disinfectants (phenols)/cold sterilants for glutaraldehyde and chlorines.

CAD/CAM Block Organizer




CAD/CAM BLOCK ORGANIZER

(50Z417) Color Code: E, J, N, Q, S, T, U

Provides a professional platform to present and discuss same-day restoration. Fits all CAD/CAM milling blocks including PrograMill®, CEREC®, PlanMill™ Blocks. Store in drawer or cabinet.

Dimensions: 14.13cm x 4.76cm x 19.68cm







E-Z Access Shelf & Accessories

E-Z ACCESS SHELF

(20Z420) Color code: A, G, I

E-Z mounting to standard equipment post $(4.44 \, \text{cm} \times 5.71 \, \text{cm})$; No tools needed. Will not damage post's finish. E-Z adjustment to multiple positions Contoured edge to contain spills and materials. Stain resistant finish; Surface disinfect only. Weight capacity of 11.34kg Rotates 360°.

Standard E-Z Access Shelf as is accommodates 4.44cm x 5.71cm post.

Workable tray area: 34.92cm x 26.03cm





E-Z ACCESS SHELF ADAPTERS

_ _

E-Z Access Shelf Adapter for Smaller Posts

(20Z422) Color code: A, G, I

Adapts to posts 3.49cm x 6.98cm

Dimensions: 5.4cm x 4.76cm x 1.9cm

E-Z Access Shelf 8.89cm Bolt Adapter Kit

(20Z437) Color code: A, G, I

Accommodates post sizes from 5.71cm x 6.03cm

Dimensions: 10.16cm x 5.08cm x 5.08cm

E-Z Access Shelf 11.43cm Bolt Adapter Kit

(20Z438) Color code: A, G, I

Accommodates post sizes from 6.03cm x 6.98cm

Dimensions: 11.43cm x 5.08cm x 5.08cm



E-Z ACCESS SHELF CUP HOLDER SNAP-ON

(20Z436I) Gray only

Fits all standard cups. No tools needed.

Dimensions: 7.94cm \times 6.35cm \times 2.54cm





E-Z ACCESS SHELF 20.32CM EXTENSION ARM

(20Z428) Color code: A, G, I

Enables shelf to extend over patient. Shelf has strength to hold 11.34kg when connected to arm. May be used with Wall Mount, Surface Mount, or directly to post.

Dimensions: 25.4cm x 10.16cm x 11.43cm





E-Z ACCESS SHELF WALL OR CABINET MOUNT

(20Z426) Color code: A, G, I

Mount on horizontal flat surface only.

Dimensions: 16.51cm x 14.60cm x 12.38cm

S



E-Z ACCESS SHELF SURFACE MOUNT

(20Z424) Color code: A, G, I

Mount on horizontal flat surface only.

Dimensions: 9.52cm x 9.52cm x 7.3cm

S

Cabinet Trays & Lab Pans

CABINET TRAYS S K





No 16A - Hand Instruments (Regular Well) (20Z201) Color code: A, B

Dimensions: 20cm x 9.52cm x 2.38cm



No 16 - Hand Instruments (Deep Well) (20Z202) Color code: A, B

Dimensions: 19.68cm x 9.52cm x 2.54cm



No 17 - Clamps, Retainers, etc. (20Z203) Color code: A, B

Dimensions: 20cm x 9.52cm x 1.9cm



No 18 - Finishing Strips (20Z204) Color code: A, B

Dimensions: 20cm x 9.52cm x 1.9cm



No 19 - Pliers, Scissors, etc. (20Z205) Color code: A, B

Dimensions: 20cm x 9.52cm x 1.9cm



No 20 - Discs, Wheels, Stones (20Z206) Color code: A, B

Dimensions: 20cm x 9.52cm x 1.9cm



LAB PANS

Organizes and transports impressions, models, dentures, and crown and bridge appliances. Stainless steel spring clip holder. Color-coded.

Lab Pan

(17Z101) Color code: A, B, D, G, M

Dimensions: 19.05cm x 14.6cm x 17.78cm

Wall-Hanging Lab Pan (17Z102) Color code: A, B, D, G, M

Dimensions: 19.68cm x 14.6cm x 17.78cm



WALL STRIP LAB PAN

(17Z104) Color code: A, B, D, G, M

Drilled and countersunk with screws and wall anchors. 45.72cm Wall Strip. Holds 3 Lab Pans.

Dimensions: 45.72cm x 2.22cm x 0.95cm

HVE & SE ACCESSORIES

Saliva Ejector Valve (with Shut-Off) 5 | 本」

- Quick disconnect with swivel
- Rubber tip with protective screen
- · Provides variable control from full on to full off







SALIVA EJECTOR SCREENS

Single use

Dimensions: 0.95cm Dia. x 0.95cm

25pk (50Z952S)

250pk (50Z952P)



Saliva Ejector Valve - Lever

SALIVA EJECTOR VALVE WITH LEVER ON/OFF CONTROL

(23E365) Gray only

Dimensions: 7.62cm x 1.59cm x 1.9cm

SALIVA EJECTOR VALVE LEVER ONLY (23E364) Gray only

Dimensions: 1.59cm x 0.48cm x 1.59cm

SALIVA EJECTOR REPLACEMENT TIP (23E363) Gray only

Dimensions: 1.59cm x 1.59cm x 2.54cm

Dispensers



COTTON ROLL DISPENSER (28R800) Color code: A, B

Heavy-duty, no-skid base. Dispenses up to three rolls at a time. Extra sturdy construction.

Dimensions: 10.64cm x 5.5cm x 10.16cm





MATRIX DISPENSER (55Z200) Color code: A, B

One roll of 0.95cm x 11m Matrix clear tape 0.002 gauge. Allows matrix to be contoured (u-shaped) as it dispenses. Stainless steel cutter allows for desired length. Disposable.

Dimensions: 5.4cm x 5.71cm x 1.59cm





PELLET DISPENSER

Spring loaded – pushes pellets to the surface.

Dispenser for Small Pellets Sizes 3, 4, 5 (28R825) White only

Dimensions: 4.76cm x 4.76cm x 5.24cm

Dispenser for Medium Pellets Sizes 0, 1, 2 (28R850) White only

Dimensions: 4.76cm x 4.76cm x 5.24cm

E-Z Storage



E-Z STORAGE TUB SYSTEM

Hinges allow cover to stay open on its own! Each tub comes with one vertical divider and one half-depth adjustment divider. Hinged cover comes in white or clear.

E-Z Storage Organizer (White Cover) (29R605) Rack sold separately

Dimensions: 14.92cm x 12.38cm x 14.6cm

E-Z Storage Organizer (Clear Cover) (29R610) Rack sold separately

Dimensions: 14.92cm x 12.38cm x 14.6cm

Replacement Hinges (Pair) (29R627)

Dimensions: 2.22cm x 1.27cm x 0.3175cm



Rack holding angled Tub

S

E-Z STORAGE SINGLE RACK (29R650)

Holds storage organizer Tub. Tub can be placed flat or at an angle for easy access of supplies. Wall or cabinet mounted with optional mounting bracket. Sturdy metal construction; plastic coated. White, non-slip feet.

Dimensions: 17.78cm x 12.7cm x 12.06cm



E-Z STORAGE MOUNTING BRACKET (29R651)

White metal bracket for the E-Z Storage Tub single rack. Pre-drilled holes for cabinet or wall mount.

Dimensions: 13.02cm x 2.54cm x 1.27cm

Handi-Hopper





HANDI-HOPPER (WITH MOUNTING BRACKET)

Small plastic waste receptacle for disposal of contaminated waste at chairside, keeping tray clean. Slotted mounting bracket (included) for easy removal and replacement. To use for infectious waste, order red (Biohazard) liners. To use for other waste, order clear liners.

Handi-Hopper (20Z404) Color code: G, I

Dimensions: 17.78cm x 10.16cm x 16.19cm

Handi-Hopper Bracket (20Z408) Color code: G, I

Dimensions: 5.71cm x 0.63cm x 2.54cm

mount this mini waste receptacle within arm's-reach of chairside for convenient disposal of contaminated waste



HANDI-HOPPER LINERS

Handi-Hopper Liners 100/Box (20Z405) Red Biohazard Single use

Handi-Hopper Liners 100/Box (20Z410) Clear Single use



WASTE-DROP RING

(22Z100) Color code: A, G, I Waste-Drop Ring for cabinet top. Fits 12.7cm hole. Made of self-lubricating, easy-to-clean resin.

Dimensions: 14.29cm Dia. x 7.62cm

S



Cuspidor/Dispos-A-Screen Traps



25PK | 144PK | DISPOS-A-SCREEN TRAPS

Avoid contact with potentially hazardous material. Improves performance of suction system. Single use.

Blue - 5.4cm 25pk | (51Z045S) 144pk | (51Z045V)

Adec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus MDT, Schein.

Green - 1.88cm 25pk | (51Z123S) 144pk | (51Z123V)

Ampco, Belmont, Biotec, Chayes, Dabi Atlante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma, and others. Gray - 6.98cm 25pk | (51Z070S) 144pk | (51Z070V)

Midwest 210 Models.

Yellow – 5.71cm 25pk | (51Z040S) 144pk | (51Z040V)

Adec, model with metal canisters.

Beige - 7.3cm 25pk | (51Z060S) 144pk | (51Z060V)

Den-Tal-EZ models with metal canister.

Mauve – 3.17cm 25pk | (51Z080S) 144pk | (51Z080V)

Adec, Cascade, and Performer.

Coral – 6.67cm 25pk | (51Z050S) 144pk | (51Z050V)

Pelton & Crane Spirit Unit with Plastic Canister.



25PK | 144PK CUSPIDOR SCREENS

Single use.

Cuspidor Screens Large Universal 6.35cm (gray) 25pk | (51Z205S) 144pk | (51Z205V)

Adec, Ampco, Belmont, Dabi Atlante, Dansereau, DCI, Den-Tal-EZ, Forest, Knight, Marus MDT, Proma, Ritter, Schein.

Cuspidor Screens Small – 4.13cm (white) 25pk | (51Z200S) 144pk | (51Z200V)

Adec, Cascade, and Performer cuspidor bowls.

Needle Capper



s |淅|

NEEDLE CAPPER

One-handed capping and recapping. Syringe is ready for use at convenient position in Needle Capper. Base of device is secured to tabletop or tray with adhesive strip (included). Needle is uncovered at the same time syringe is withdrawn from Capper. After use, syringe is replaced for recapping. Needle Capper available in Standard and Long sizes. Either size tube will fit standard base.

Standard Needle Capper (27R700) White only

Includes: base, one tube, adhesive tape.

Dimensions: 5.56cm x 3.17cm x 3.17cm

Long Needle Capper (27R704) White only

Includes: base, one tube, adhesive tape.

Dimensions: 5.63cm x 3.17cm x 4.44cm

Mighty Mixer Bowls, Spatulas, and Organizers



MIXING STICKS 100PK (50Z521) White only

Single use.

Dimensions: 11.43cm long

DELRIN SPATULA 10PK (50Z523)

Flexible flat blade. Non-stick. Resistant to chemicals. Great for mixing cements, acrylics, etc.

Dimensions: 11.75cm x 0.48cm









CEMENT SPATULA (50Z502)

White, plastic, double-end spatula. For fillings, cement, and pastes. Selflubricating plastic for easy cleaning.

Dimensions: 15.19cm x 0.63cm x 0.95cm







MIXING WELL 144PK (20Z496V) White only

Single use.

Dimensions: 11.43cm long







MIGHTY MIXER SPATULAS

Self-lubricating for easy cleaning. Broad, medium flex blade. Ideal for alginates and plastic.

Dimensions: 18.57cm x 0.48cm x 3.17cm





Single Spatula (50Z501) White only



Assorted 4 Pack (50Z501ASTD) 1 of each color: J, N, P, and R

Dimensions: 18.57cm x 0.16cm x 3.17cm



MIGHTY MIXER BOWLS 36PK (50Z505) White (50Z505ASTD) Assorted

Disposable, eliminating cleanup time. Fits securely into mixing bowl without slippage. Saves wear on mixing bowl. Free-standing for use with or without mixing bowl. Capable for use with electric Alginate mixers. Single use.

Dimensions: 12.38cm Dia. x 7.94cm



IMPLANT ORGANIZER

(20Z490B) Blue only

Quick and easy recognition. Individual compartments for implants. Shape of the organizer matches the arch of the mouth for implant placement. Separate compartments prevent implant appliances from mixing together. Lubricated plastic enables tooth number to be written in the designated area. Cover helps keep materials secure.

Dimensions: 16.19cm x 0.63cm x 0.95cm





RESTORATIVE ORGANIZER

(20Z495R) Neon purple only

Compartments shaped like the arch of the mouth for veneers, inlays, composite, and onlays. Disposable Mixing Wells included for bonding agent and cement. Accessory area for brushes and materials. Lubricated plastic enables tooth number to be written in the designated area. Light protective cover can be used as a light shield and to prevent materials from setting early. Includes 25pk of Mixing Wells.

Dimensions: 17.46cm x 11.75cm x 1.59cm

S

RELIABILITY CUSTOMER FOCUS PASSION INTEGRITY CREATIVITY FIND SOLUTIONS MADE IN THE USA





sales@zirc.com

20Z500C-EU REV B 0123