

PROMOTIONAL OFFERS*

Valid from 3th of February till 30th of June 2025

Only for participating dealers Benelux.

All prices are in €, exclusive VAT.

newtron p5XS

1

Ref. F00246



P5XS b-LED incl. with the following tips:
Ref. F62100

~~€2788,25~~

€1671

1S

Ref. F00245



H3

Ref. F00369



N°4

Ø22


Ref. MP3622PH

NOW INCLUDED

PURE REFLECT®
with a value
of €54,89

Ref. MP3622PH

Newtron® handpiece



ref. F12281 - Autoclavable

~~€600,77~~

€245

Newtron® LED handpiece



ref. F12609

~~€831,57~~

€385

Newtron® Slim B. LED handpiece



- with blue LED ring: ref. F12900
- with white LED ring: ref. F12905

~~€831,57~~

€467

Hygiene kit

ref. F00935



~~€516,23~~

€370

Perio kit

ref. F00936



~~€654,60~~

€469

Implant Protect kit

ref. F02120



~~€698,16~~

€533

Endo Success Retreatment kit

ref. F00737



~~€667,28~~

€494

piezotome cube

Dynamic power adaptation for
safe bone surgeries

Ref. F50100

Delivered with a multifunction footswitch,
a Piezotome® CUBE LED handpiece,
1 handpiece holder, an Essential II kit
(6 tips BS1S, BS4, SL1, SL2, SL3 et LC2),
a Piezotome® dynamometric wrench,
1 irrigation pole, 2 autoclavable irrigation
lines, 5 clips and 15 perforators.

~~€7.423,42~~

€4499



Intraoral surgery ultrasonic control console.
Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

implant center 2

3 modes:

- Piezotome® (bone surgery).
- Newtron® (perio, endo,...).
- I-Surge™ (implant motor).

Ref. F27200 - Delivered with a LED motor, a LED motor cord, a progressive footswitch, 2 brackets, 2 handpiece holders, 2 sterilizable irrigation lines and 30 perforators.

~~€14.183,77~~

€7643



The alliance of technologies
for powerful and secure
bone surgeries

Ultrasonic control console for conventional dentistry and intraoral dentistry and motor control console for dental implantology
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France



Bone Surgery™ II Kit*

Ref. F87509 - Tips BS1S II, BS2L II, BS2R II, BS4 II, BS5 II, BS6 II.

~~€1359,48~~

€868



Sinus Lift II Kit*

Ref. F87519 - Tips SL1 II, SL2 II, SL3 II, SL4 II, SL5 II.

~~€1051,63~~

€788



Intralift™ II Kit*

Ref. F87536 - Tips TKW1 II, TKW2 II, TKW3 II, TKW4 II, TKW5 II (Ø 3 mm).

~~€1334,35~~

€868



Extraction™ II Kit*

Ref. F87546 - Tips LC1 II, LC2 II x2, LC2L II, LC2R II, Ninja® II.

~~€1032,24~~

€788

* Kits are delivered with an autoclavable metal tip holder and a universal wrench.
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

Accessoires

Promotion valable only for ImplantCenter2

Newtron® LED handpiece



Ref. F12609
For Piezotome® 2 and ImplantCenter™ 2

Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

Newtron® LED cord



Ref. F12830
For Piezotome® 2 and ImplantCenter™ 2

~~€1.811~~
€1186

Universal irrigation lines Sterile x10



€115
as from 4 pieces

Ref. F59905 - Single-use with integrated cassette and flow control system.

Class Is medical device (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

PIEZOTOME® LED handpiece + cord



Ref. F12815

Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

~~€1.592,18~~
€1034



Universal II Kit*

Ref. F87539 - Tips BS6, LC1, LC2, CE1, CE3 et PZ1.

~~€1359,48~~
€868



Essential II Kit*

Ref. F87528 - Tips BS1S II, BS4 II, SL1 II, SL2 II, SL3 II, LC2 II.

~~€1359,48~~
€868



Crest Splitting II Kit*

Ref. F87567 - Tips CS1 II, CS2 II, CS3 II, CS4 II, CS5 II, CS6 II.

~~€880,10~~
€652



Crown Extension II Kit*

Ref. F87554 - Tips BS6 II, CE1 II, CE2 II, CE3 II.

~~€684,40~~
€493

* Kits are delivered with an autoclavable metal tip holder and a universal wrench.

servotome

Delivered with a single footswitch,
a box with 10 electrodes,
a conductive bracelet.

Ref. F57210

~~€1.976,45~~

€1353



Electrodes

Fulguration & coagulation




Coagulating incision









Excision



Incision



| |  |  |
|---|--|---|
|  | 6 - 8 | 3 - 4 |
|  | 6 - 8 | 3 - 4 |
|  | 6 - 8 | 3 - 4 |
|  | 5 - 7 | 1 - 3 |

air-n-go easy

All-in-one polisher for easy procedures

- AIR-N-GO® easy with KaVo* adapter Ref. F10100
- AIR-N-GO® easy with Sirona* adapter Ref. F10101
- AIR-N-GO® easy with W&H* adapter Ref. F10102
- AIR-N-GO® easy with Bien-Air* adapter Ref. F10103
- AIR-N-GO® easy with Midwest* adapter Ref. F10105



with a convertible AIR-N-GO® easy handpiece, a Supra nozzle, 1 CLASSIC powder bottle neutral (250 g), 1 AIR-N-GO® easy maintenance kit (replacement o'rings, syringe of silicone grease with brushes, cannula, cleaning drill).

*KaVo, Sirona, W&H, Bien-Air and Midwest are registered trademarks belonging respectively to the company of the same name.

air-n-go perio



air-n-go easy + air-n-go perio > ~~€2.070,30~~

€1258

Kit PERIO easy

with a AIR-N-GO® easy nozzle and a bottle of AIR-N-GO® PERIO glycine powder 160g of AIR-N-GO® PERIO glycine powder.

Ref. F10120

Dental air-polishing handpiece - Class IIa medical device (GMED) - CE 0459.

Manufacturer: SATELEC® a company of ACTEON® group - France

AIR-N-GO® "CLASSIC" powder

4 bottles 250 g

- Peppermint Ref. F10250
- Lemon Ref. F10251
- Raspberry Ref. F10252
- Cola Ref. F10253
- Neutral Ref. F10254
- Assortment Ref. F10255

€79

each

€59

as from
4 pieces



air-n-go perio



AIR-N-GO® "PERIO" powder

Ref. 10270

3 bottles 160 g

~~€126,97~~

€91

AIR-N-GO® CLASSIC: Class IIa medical device - (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

AIR-N-GO® PERIO: Powder for subgingival treatment of periodontitis and peri-implantitis - Class IIa medical device - (GMED) - CE 0459

Manufacturer: SATELEC® a company of ACTEON® group - France

miniled standard

~~€1.095,56~~

€708

Ref. F02507

- **Fast**

Only 6 to 12 sec to polymerize
2mm composite.

- **Power**

1250mW/cm²
(Ø7,5mm light guide).

2000mW/cm²
(Ø5,5mm light guide - optional).

- **Safety**


Base-charger includes a radiometer for light power
intensity check-up.

- **Optimal result**

3 modes to fit all clinical needs.

Highly Adaptable


The alliance of power & safety



FAST MODE

Fast and easy


Full power during 10 sec
(bip at 5 sec) on regular
depth.



PULSE MODE

Ally rapidity & safety

10 flashes of 1 sec
(bip after 5 flashes)
with times off to
preserve dental pulp.



RAMPING MODE

No stress cracking

Progressive power
increase during 20 sec
(bip every 5 sec) for
a complete and safe
curing.

Simply Powerful

Outperforms in
power, efficiency and speed

~~€1.407,36~~

€868

Ref. F05214



miniled supercharged

- **Visibility**

Luminous screen with customizable modes and times, count curing time and load indicator.



- **Power**

2000mW/cm²
(Ø7,5mm light guide).

3000mW/cm²
(Ø5,5mm light guide - optional).

- **Optimal result**

3 modes to fit all clinical needs.



FAST MODE

The choice of rapidity


Full power during 3, 4, 5 or 10 sec, according to the composite thickness.



PULSE MODE

Limits temperature rise

5 or 10 flashes of 1 sec for a polymerization in security.



RAMPING MODE

Avoids composites retraction

Progressive power increase during 9 or 15 sec to ensure perfect sealing of restoration.

miniled ortho 2

~~€1.457,10~~

€868

Ref. F05220

- **Customizable signal**

An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec.

- **Power**

3000mW/cm²
(Ø5,5mm light guide).

- **Long lasting battery**

MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day.

- **Ortho friendly**

Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch.

- **Efficacy & safety**

Cure any type of material used orthodontics efficiently and quickly.

No noisy and unhygienic fan needed to cool the lamp.

Specialized Orthodontic Treatments

Cure in a simple gesture, in few seconds from
few brackets to a whole quadrant



etik collagene

Hemostatic surgical dressing

Etik™ Collagene

25 unit-packaged sterile cylinders (in blister packs)

Ref. 281 970

~~€90,45~~

€63



Class III medical device (GMED) - CE 0459 - Manufacturer: PRODUITS DENTAIRES PIERRE ROLLAND - ACTEON® Group - France

softly™

Pre-measured
capsules amalgamator

Ref. 6.01.00112

~~€496,85~~

€330



Medical device classe I - CE - Manufacturer: SATELEC® Groupe ACTEON® (France)

Discover our RISKONTROL® range

Adapt your tips to the taste

of your patients



art
by **riskontrol**

~~€35,20~~

€26



Ref. 201 700 - Mint



Ref. 201 701 - Mandarine



Ref. 201 702 - Liquorice



Ref. 201 703 - Anise



Ref. 201 704 - Black currant

Non-sterile disposable tips (x 250)

When you buy
10 boxes 1 adaptor FOC

Add color

to your patients' care experience



~~€35,20~~

€26

classic
by **riskontrol**



Ref. 201 664 - White



Ref. 201 665 - Blue



Ref. 201 666 - Yellow



Ref. 201 667 - Green



Ref. 201 668 - Pink

Non-sterile disposable tips (x 250)

Protect entirely your air/water syringe

thanks to the plastic sheath



total protec
by **riskontrol**

~~€31,90~~

€24



Ref. 201 750



Ref. 201 751

Non-sterile disposable tips with plastic sheath (x 100)

Choose a totally sterile

workflow



~~€33,00~~

€24



sterile
by **riskontrol**




Ref. 201 770

Sterile disposable tips (x 96)

Classic by RISKONTROL®, Art by RISKONTROL®, Total Protec by RISKONTROL®: Class I medical device. Sterile by RISKONTROL®: Class Is medical device. CE0459 (GMED). Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of these devices is in the Product Feature Summary available on the laboratory's website. Manufacturer: PRODUITS DENTAIRE PIERRE ROLLAND - ACTEON® Group - FRANCE.



A TRIPLE ACTION



Kaolin clay



Aluminum chloride

1 ASTRINGENT EFFECT

helps to completely dry the sulcus with no blood or fluids.

2 OPEN THE SULCUS

effectively with no pain or damage to tissues thanks to the viscosity carefully calculated with our clinical expertise.

3 PRESERVE EPITHELIAL ATTACHMENT*

for a reduced risk of gingival recession and bone resorption.

**Canula must not be inserted into the sulcus*

**expasyl
exact**



~~€140,81~~

€75

Paste for temporary sulcular opening
20 capsules (0,3g) Expasyl™ Exact.



expasyl

~~€238,04~~

€99

Paste for temporary sulcular opening
20 capsules (1g) Expasyl™.



~~€106,36~~

€60

Straight - Ref. 261005
Bended - Ref. 261015